
	TITLE: Approved Support Surface Cleaners		DOCUMENT #: REF000115
	REVISION: A	REVISION DATE: 10 July 2023	DCR #: DCR-01453
	OWNING DEPARTMENT: Engineering		PAGE: 1 of 9


Printed copies must always be compared against the electronically filed controlled document to ensure current version.

Manufacturer	Product Description	Concentration
(Various)	Sodium Hypochlorite	1.3% (1.3:100 water)
(Various)	Hydrogen Peroxide 5%	Ready to use
(Various)	Isopropyl Alcohol 70%	Ready to use
(Various)	Trix (Soap)	Ready to use
3M	TB Quat Disinfectant	Ready to use
ABC Compounding Co., Inc.	Aero Stain Away	Follow chemical manufacturer's specifications
Agar Cleaning Systems Pty. Ltd.	Chloradet	0.6% Chloradet (1:160 water)
Alco Mandala Products	DC-7	Follow chemical manufacturer's specifications
Alco Mandala Products	Sanni-Rinse	Follow chemical manufacturer's specifications
Ameriderm Laboratories Limited	Care Source Disinfectant	Follow chemical manufacturer's specifications
Anios Laboratories	Surfa 'Safe	Follow chemical manufacturer's specifications
Anti-Slip Anywhere	Multi Kleen	Follow chemical manufacturer's specifications
Apollo Corporation	CID A-L II Quaternary Disinfectant	Follow chemical manufacturer's specifications
Apollo Corporation	Hygena	Follow chemical manufacturer's specifications

	TITLE: Approved Support Surface Cleaners		DOCUMENT #: REF000115
	REVISION: A	REVISION DATE: 10 July 2023	DCR #: DCR-01453
	OWNING DEPARTMENT: Engineering		PAGE: 2 of 9


Printed copies must always be compared against the electronically filed controlled document to ensure current version.

Manufacturer	Product Description	Concentration
Apollo Corporation	Therasol	Follow chemical manufacturer's specifications
Apollo Corporation	Turbo Clean	Follow chemical manufacturer's specifications
ArjoHuntleigh	Disinfectant Cleanser IV	Follow chemical manufacturer's specifications
Bioscience Fresche	Surface Plus	Follow chemical manufacturer's specifications
Brighton Professional	Hepastat 256	Follow chemical manufacturer's specifications
Buckeye International, Inc.	Eco E23 S23	256:1 dilution with water
Buckeye International, Inc.	Sanicare TBX	Follow chemical manufacturer's specifications
Buckeye International, Inc.	Terminator One Step Disinfectant	Follow chemical manufacturer's specifications
Buckeye International, Inc.	Sanicare Disinfectant Wipes	Ready to use
Certo	Non Acid Bathroom Cleaner	Follow chemical manufacturer's specifications
Clinell	Detergent Wipes	Follow chemical manufacturer's specifications
Clinell	Hand & Surface Sanitizer	Follow chemical manufacturer's specifications
Clinell	Universal Disinfectant Wipes	Follow chemical manufacturer's specifications
Clinell	Sporicidal Wipes	Ready to use
Clinell	Universal Wipes	Ready to use

	TITLE: Approved Support Surface Cleaners		DOCUMENT #: REF000115
	REVISION: A	REVISION DATE: 10 July 2023	DCR #: DCR-01453
	OWNING DEPARTMENT: Engineering		PAGE: 3 of 9


Printed copies must always be compared against the electronically filed controlled document to ensure current version.

Manufacturer	Product Description	Concentration
Clorox	Bleach Wipes (Sodium Hypochlorite 0.1 - 1%)	Follow chemical manufacturer's specifications
Clorox	Dispatch Disinfectant	Follow chemical manufacturer's specifications
Clorox Disinfecting Wipes	Clorox Disinfecting Wipes	Ready to use
Clorox HEALTHCARE	Fuzion	Follow chemical manufacturer's specifications
Clorox HEALTHCARE	Hydrogen Peroxide Cleaner Disinfectant Wipes	Ready to use
Delaval	Delaval Acid Wash	5% Delaval (1:20 water)
Dupont	Virkon	5% Virkon (1:20 water)
Ecolab	Chlorwhite	1% Chlorwhite (1:100 water)
Ecolab	A-456-II	1.6% A-456-II (1:64 water)
Ecolab	Actichlor Plus	2000PPM & 10,000PPM
Ecolab	Oasis Pro 12 Neutral Cleaner	5% Oasis Pro 12 (1:20 water)
Ecolab	Spring Clean	5% Spring Clean (1:20 water)
Ecolab	Asepti Active Liquid	Follow chemical manufacturer's specifications
Ecolab	Incidin Extra	Follow chemical manufacturer's specifications

	TITLE: Approved Support Surface Cleaners		DOCUMENT #: REF000115
	REVISION: A	REVISION DATE: 10 July 2023	DCR #: DCR-01453
	OWNING DEPARTMENT: Engineering		PAGE: 4 of 9


Printed copies must always be compared against the electronically filed controlled document to ensure current version.

Manufacturer	Product Description	Concentration
Ecolab	Incidin Foam	Follow chemical manufacturer's specifications
Ecolab	Incidin Plus	Follow chemical manufacturer's specifications
Ecolab	Oxycide	Follow chemical manufacturer's specifications
Ecolab	Asepti Active Wipes	Ready to use
Ecolab Cleantec	Nature Clean	1:38 water
Hartmann	Sterillium	Follow chemical manufacturer's specifications
Helix Solutions	Chlor-Clean	1000PPM (1 tablet per 1 liter water)
Hu-Friedy Mfg. Co., LLC.	Sporeclear RTU Spray or Wipes	Diluted at 8% into water
Hunters Products	TRI-DET	1:80 water
Kimberly-Clark	Isowipes	Ready to use
Legge Systems	Legcide HB Plus	Follow chemical manufacturer's specifications
Livingstone International (PTY) LTD	LIV-Wipes Jumbo	Ready to use
MEDCHEM	BAK 1:750	Follow chemical manufacturer's specifications
Medizar	Multi-Surface Sanitiser Wipes	Ready to use

	TITLE: Approved Support Surface Cleaners		DOCUMENT #: REF000115
	REVISION: A	REVISION DATE: 10 July 2023	DCR #: DCR-01453
	OWNING DEPARTMENT: Engineering		PAGE: 5 of 9


Printed copies must always be compared against the electronically filed controlled document to ensure current version.

Manufacturer	Product Description	Concentration
Medline Industries Inc.	Micro-Kill Q3 Concentrated Disinfectant, Cleaner & Deodorizer	Follow chemical manufacturer's specifications
Medline Industries Inc.	MICRO-KILL AF2 Wipes	Ready to use
Medline Industries Inc.	MICRO-KILL Bleach Germicidal Wipes	Ready to use
Medline Industries Inc.	MICRO-KILL Bleach-Free, Alcohol-Free Wipes	Ready to use
Medline Industries Inc.	MICRO-KILL One Wipes	Ready to use
Medline Industries Inc.	Touch Screen Cleaning Wipes	Ready to use
Metrex Research	CaviCide	Follow chemical manufacturer's specifications
Metrex Research	CaviCide1	Follow chemical manufacturer's specifications
Metrex Research	CaviWipes1	Ready to use
Metrex Research	CaviWipes	Ready to use
Milton	Anti-Bacterial Tablets	1 Tablet:2L water
Noble Pine Products Co.	Steri-Fab	Follow chemical manufacturer's specifications
OdorXit Naturals	Odor Eliminator Spray	Ready to use
Pal International	PAL Detergent Wipes	Ready to use
Pal International	PAL Disinfectant Wipes	Ready to use

	TITLE: Approved Support Surface Cleaners		DOCUMENT #: REF000115
	REVISION: A	REVISION DATE: 10 July 2023	DCR #: DCR-01453
	OWNING DEPARTMENT: Engineering		PAGE: 6 of 9


Printed copies must always be compared against the electronically filed controlled document to ensure current version.

Manufacturer	Product Description	Concentration
PCS Process Cleaning Solutions	PCS 7000	Follow chemical manufacturer's specifications
PCS Process Cleaning Solutions	PCS Spray & Wipe Hand and Surface Cleaner	Ready to use
PDI	Super Sani Cloth Germicidal Disposable Wipes	Ready to use
PreVasive	BAC Botanical Antimicrobial Cleaner	Follow chemical manufacturer's specifications
PreVasive	Noroxycdiff	Follow chemical manufacturer's specifications
Pro Restore Products	Microban X-590 Institutional Spray	Follow chemical manufacturer's specifications
Pure Green, LLC	Puregreen 24 (Pure Green)	Follow chemical manufacturer's specifications
Quantum Pacific Limited	Actimax	Follow chemical manufacturer's specifications
Reynard Health Supplies	Reynard Neutral Detergent Wipes	Ready to use
Reynard Health Supplies	Reynard Surface Sanitizing Wipes	Ready to use
Reynard Health Supplies	Surface Sanitising Wipes	Ready to use
sBioMed, LLC	Steriplex SD Disinfectant Sporicide Wipes (Without Bleach)	Follow chemical manufacturer's specifications
Schulke & Mayr GMBH	Kodan Tinkture Forte	Follow chemical manufacturer's specifications
Schülke & Mayr GmbH	Ocetenisept	Follow chemical manufacturer's specifications

	TITLE: Approved Support Surface Cleaners		DOCUMENT #: REF000115
	REVISION: A	REVISION DATE: 10 July 2023	DCR #: DCR-01453
	OWNING DEPARTMENT: Engineering		PAGE: 7 of 9


Printed copies must always be compared against the electronically filed controlled document to ensure current version.

Manufacturer	Product Description	Concentration
SCIENCECompany	Ammonium Chloride	Follow chemical manufacturer's specifications
Sealed Air Diversey Care	Virex II 256	1:256 water
Sealed Air Diversey Care	Stride Citrus Smartdose	1:750 water
Sealed Air Diversey Care	Divercleanse	3.1% Divercleanse (3.1:100 water)
Sealed Air Diversey Care	Room Care R2L	Follow chemical manufacturer's specifications
Sealed Air Diversey Care	Stride HC (JF)	Follow chemical manufacturer's specifications
Sealed Air Diversey Care	Suma San Conc D4A	Follow chemical manufacturer's specifications
Sealed Air Diversey Care	Suma San Conc D4A	Follow chemical manufacturer's specifications
Sealed Air Diversey Care	Virex Plus	Follow chemical manufacturer's specifications
Sealed Air Diversey Care	Oxivir TB Liquid and Wipes	Ready to use
Spartan Chemical Company Inc.	Super HDQ Neutral	256:1 dilution with water
Spartan Chemical Company Inc.	BNC-15	Follow chemical manufacturer's specifications
Spartan Chemical Company Inc.	GS Neutral Disinfectant #3502	Follow chemical manufacturer's specifications
Spartan Chemical Company Inc.	hdq C2	Follow chemical manufacturer's specifications

	TITLE: Approved Support Surface Cleaners		DOCUMENT #: REF000115
	REVISION: A	REVISION DATE: 10 July 2023	DCR #: DCR-01453
	OWNING DEPARTMENT: Engineering		PAGE: 8 of 9

Printed copies must always be compared against the electronically filed controlled document to ensure current version.

Manufacturer	Product Description	Concentration
STERI-7	S-7 Xtra Wipes	Ready to use
Sterri-Matt	Sterri-Matt (Pure) Foam	Ready to use
Sulco LTD	Mediwipes	Ready to use
Synergy Health	Azomax Multi Surface Wipes	Ready to use
Synergy Health	Azowipe Hard Surface Wipes	Ready to use
Texwipe	TX 3270	Ready to use
Tristel Solutions LTD	Distel Surface Disinfectant (formerly Trigene)	Follow chemical manufacturer's specifications
True Blue Chemicals	Use All	Follow chemical manufacturer's specifications
Velocity Chemicals	Eco-Pine	5% Ecopine (1:20 water)
Verna Care	Tuffie 5 Cleaning Wipes	Ready to use
Wexford Labs	Wexcide 128	Minimum 128:1 dilution with water
Wexford Labs, Inc.	CleanCide Wipes	Ready to use
Whiteley Medical	Lemex	Follow chemical manufacturer's specifications
Whiteley Medical	Odacon Spray	Follow chemical manufacturer's specifications
Whiteley Medical	Viraclean	Follow chemical manufacturer's specifications

	TITLE: Approved Support Surface Cleaners		DOCUMENT #: REF000115
	REVISION: A	REVISION DATE: 10 July 2023	DCR #: DCR-01453
	OWNING DEPARTMENT: Engineering		PAGE: 9 of 9

Printed copies must always be compared against the electronically filed controlled document to ensure current version.

Manufacturer	Product Description	Concentration
Wipesplus	Hospital/Surgical Wipes	Ready to use

NOTES:

To avoid premature degradation of the support surface fabrics, never use an abrasive cloth, sponge or other similar devices on the surface. **Wipe off cleaning product residue with a clean, soft cloth and rinse with fresh water immediately after use.** Allow the surface to dry completely before use. Failure to follow these procedures may reduce the lifespan of the fabric.

Cleaners/disinfectants can cause discoloration/bleaching on some fabrics. To avoid discoloration, always check on a small area first. Information contained here is purely related to the compatibility of cleaners with Agiliti support surface materials and in no way relates to the efficacy of any cleaner for infection control or the safety of any cleaner for skin contact.

All approvals are based on strictly following the chemical manufacturer's recommended cleaning procedures, including dilution ratios (if noted).

If a particular cleaner/disinfectant is not listed, contact Agiliti for further recommendation.