


The 4 Technologies Behind NPT4™



Top Cover with AirSpace™ 3D Mesh:

Sizewise is known for its fabrics, and this three-ply top cover is like nothing we've ever built. The extreme polycarbonate fabric, AirSpace™ 3D Mesh, and Mattress Asset Protection (MAP) Liner™ work together to keep the patient dry and reduce shear.



MAP Liner™:

Hospital cleaning protocols are no match for NPT4's high-tech design and construction. Our exclusive MAP Liner™ not only bolsters the top cover's performance, it also provides a level of infection control and prevention you won't find elsewhere—extending the life of the mattress.



CuPro™ Copper-Infused Foam:

NPT4 introduces CuPro™ copper-infused, thermoresponsive foam to our lineup. It breathes more freely, adding to NPT4's microclimate management performance, and the copper particles help dissipate heat more efficiently, keeping the patient cool while reacting to temperature and weight.



Microclimate Management Pump:

NPT4's easy-to-use pump offers 1-in-2 alternation therapy and a layer of air from head-to-heel—providing advanced skin protection to reduce pressure injuries and ensure patient comfort.

Learn more: Call 800-814-9389
or contact your Sizewise rep