

# Evolution<sup>®</sup> Expandable

Fall Prevention Bed



*Shown with optional defined perimeter foam support surface*

Evolution<sup>®</sup> Expandable boasts an industry-leading 9" deck low and a 30" deck high for an ergonomic working environment. The bed transports patients in any high or low position, and our Rest Secure System<sup>™</sup> (RSS) comes standard. Evolution Expandable accommodates mattresses from 39-48". Safe Working Load: 850 lbs.

- *Focus on fall prevention: 9" deck low*
- *Rest Secure System<sup>™</sup> (RSS) includes bed exit alarm, protocol timer, and scale*
- *Transports in any high/low position*
- *Accommodates mattresses from 39-48"*

sizewise.com 800-814-9389

  
**SizeWise<sup>®</sup>**



- 1 — Evolution® Expandable has a deck low of 9", the lowest full-feature bed frame. Also available, Low Boy(R) Expandable, with an industry leading 7.5" low deck
- 2 — Floor Safety Mat included with rental
- 3 — Rest Secure System™ includes bed exit alarm, protocol timer, and scale

## Key Specs

### BED

Low Position – Deck to Floor	9"
High Position – Deck to Floor	30.75"
Deck Width	39", 42", or 48"
Deck Length	82"
Safe Working Load	850 lbs.

### SIDE RAILS

Type	Fold Down
Included	2 Standard; Foot Optional

### HEAD/FOOT BOARDS

Material	High-Impact Plastic
----------	---------------------

### CASTER AND BRAKING SYSTEM

Caster Diameter	4"
Dual Brake/Steer System	4"

### BED ANGLES

Head Section	0-58°
Knee Section	0-49°
Trend	15°
Reverse Trend	-15°

### OPTIONS AND ACCESSORIES

Battery Backup	Standard
CPR Release	Standard
Floor Safety Mat	Standard
Footboard Lockout Controls	Standard
Rest Secure System™ (RSS)	Standard
<i>Bed Exit Alarm, Protocol Timer, Built-in Scale</i>	
IV Pole	Optional
O2 Holder	Optional
Pillow Speaker	Optional
Seizure Pad Covers	Optional
Trapeze	Optional
6-Port Aux Outlet	Optional
Color: Mist	Standard
Color: Cocoa	Optional

### MATTRESS OPTIONS

Numerous non-powered and low air loss therapeutic surfaces available.