



BARIATRIC MOBILITY

BARIATRIC WALKER

Features

- Folding bariatric walker with increased width and depth to accommodate larger users
- Dual-release folding mechanism enables user to fold walker sides independently
- Comfortable, easy to hold handgrips
- Easy push buttons to facilitate height adjustment best suited for user
- 750 lb maximum weight capacity - EVENLY DISTRIBUTED on non-wheeled model
- 650 lb maximum weight capacity - EVENLY DISTRIBUTED when using bariatric wheel kit
- Designed to meet the requirements of HCPCS code: E0148
- Three (3) year limited warranty



WIB | Kit Includes:


- 5176 1 EA
- 5176-1 - Wheel Kit

Specifications

5176	Width between Handgrips		overall Base Width		overall Base Depth		Minimum Height		Maximum Height		Recommended User Height		Product Weight		Weight Capacity EVENLY DISTRIBUTED		Color
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	lb	Kg	
5176 non-wheeled															750	390	Black
5176 with wheels	20	50.8	23	58.4	21	53.3	32	81.3	39	99.1	5'2"—6'4"	157.5—193.0	13.5	6.1	650	294	

WARNING: Cancer and Reproductive Harm - www.P65Warnings.ca.gov

Color variations may occur due to monitor or printer settings.



+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website. www.grahamfield.com