



umano medical

MEDICAL SURGICAL BED WITH SCALE AND BED EXIT DETECTION SYSTEM



10" / 25.4 cm low height with built-in bed exit detection system to minimize falls

Seamless design to maximize infection control

Easy to use and fully featured

LOW HEIGHT • EASY TO CLEAN • EASY TO USE FULLY FEATURED




Seamless smooth surfaces maximizing infection control



Conveniently located USB power outlets to plug in mobile devices



Open architecture to simplify cleaning and mattress installation



Under bed clearance for lift and overbed table compatibility



UmanoMedical

MEDICAL SURGICAL BED WITH SCALE AND BED EXIT DETECTION SYSTEM

BED FEATURES

- 4 DC actuators
Battery back-up
4 casters with central brake
One hand release siderails
Open architecture
One step Trendelenburg
One step vascular foot position
Integrated pump holder
Transport handles on headboard and footboard
4 corner roller bumpers
Nurse controls integrated into siderails and footboard
Patient controls integrated into siderails
Patient control lockout
Boostless™ Backrest system
Backrest angle indicator
Compliant with standard IEC 60601-2-52

OPTIONAL FEATURES

- 10" / 25.4 cm low height
Built-in-bed scale system
Built-in-bed exit detection system - 2 or 3 levels of sensitivity
Remote bed exit alarm connectivity
Bed exit side view lights
Smart screen on footboard control
Trendelenburg angle displayed on footboard
Auto-contour
Programmable lowest height
Patient surface length (80" or 84")
3 position integrated bed extender : 80" / 203.2 cm 84" / 213.4 cm 90" / 228.6 cm
5" / 12.7 cm single or dual casters (optional antistatic casters)
6" dual casters
120-volt / 5A or 220-volt / 2A insulated outlets with USB mobile device power outlets
Foot end back-up control on the frame
CPR release handle
Roller type wall protection bumpers
Higher siderails (for 8" mattress)
Foot control system
Downward motion interrupt
Communication port (for upgradable bed software)
Smart Nightlights™
Bed status monitoring
One-touch button groups
30-degree head of bed
One-button chair position
Flat bed position
Backlit patient controls
Room light and/or reading light control
Nurse call

Accessories*

- Patient helper system
Upright oxygen bottle holder
IV poles (permanent & removable)
6" or 10" / 15 or 25.4 cm bed extender cushion
* To see all accessories available, visit our website.

Specifications



Table with 2 columns: Specification and Value. Rows include Product model (FL36, FL36-SC), Overall length (88", 92"-98"), Overall width (39.25", 40.5"), Weight capacity (600 lb), Height range (10"-31.5"), Patient surface positioning (Backrest, Knee gatch, etc.), Patient surface dimensions (35" x 80", etc.), Caster diameter (5"), and Bed weight (375 to 438 lb).

Electronics

Table with 2 columns: Specification and Value. Rows include Compliant with standard (IEC 60601-1 3rd edition), Volts (120 V AC, 220-240 V AC), Ampere rating (4.25A, 2.15A), Current leakage (< 300 micro amps), Frequency (60 Hz / 50 Hz), Hospital grade plug (Yes), Battery type (Sealed rechargeable lead acid 1.2Ah, 2 * 12V), Scale system accuracy (0.5%), and Scale accuracy (EC) (OIML R76-1 Class III MAX 272.5 kg).

Umano Medical reserves the right to change specifications without notice.



Umano Medical Inc.

230, boul. Nilus-Leclerc, L'Islet QC G0R 2C0 CANADA
Toll free number : 1 844-409-4030
information@umanomedical.com
umanomedical.com