

QUICK-SHIP FOAM SUPPORT SURFACES

Agiliti quick-ship options provide facilities a fast, cost-effective way to replace damaged foam support surfaces. Products can be delivered in limited quantities within five business days. Contact your Agiliti representative for production timeframe on larger quantities.



NP Trio™



NP Adjust™



NPT4™

Product Description	3-layer, 3-zone non-powered foam support surface	Vented, self-adjusting foam support surface	Vented, self-adjusting foam support surface with optional pump
Standard Size	35"W x 82"L	35"W x 82"L	35"W x 82"L
Max Patient Weight	500 lb.	500 lb.	500 lb.
Warranty: Top/Bottom Cover*	4 Years	4 Years	2 Years
Warranty: RF-Welded Seams*	4 Years	4 Years	5 Years
Product Webpage	agilitihealth.com/product/np-trio	agilitihealth.com/product/np-adjust	agilitihealth.com/product/npt4

These Agiliti-manufactured support surfaces feature construction and design details that help serve patients in a busy hospital environment:



4-way stretch top cover fabric is durable, breathable, waterproof and highly chemically resistant



CoreShield™ liner provides an additional layer of waterproof protection to the foam core (MAP Liner™ on NPT4)



RF-welded seams protect against fluid ingress, extending the useful life of the product



CuPro™ copper- and gel-infused foam helps dissipate heat from the patient/surface interface

Learn more:

Contact your Agiliti sales representative or call 800-814-9389